

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse anti Human CD2, IgG2a, Clone: [T6.3], Unconjugated, Monoclonal NMB-0021

Artikelname	Mouse anti Human CD2, IgG2a, Clone: [T6.3], Unconjugated, Monoclonal
Artikelnummer	NMB-0021
Hersteller Artikelnummer	0021
Alternativnummer	NMB-0021
Hersteller	NordicMubio
Wirt	Mouse
Kategorie	Antikörper
Applikation	FC
Spezies Reaktivität	Human
Konjugation	Unconjugated
Produktbeschreibung	Mouse anti human CD2 can be used for identification of human T cells and a subset of NK cells associated with the receptor for sheep erythrocytes rosettes expressing the 45-50,000 M.W. surface antigen....
Klonalität	Monoclonal
Konzentration	See vial for concentration
Klon-Bezeichnung	[T6.3]
Isotyp	IgG2a

UniProt	P06729
Puffer	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Quelle	Derived from hybridization of mouse Sp2/0 myeloma cells with spleen cells from BALB/c mice immunized with T Lymphocytes activated by a mixed lymphocyte culture.
Reinheit	Protein A/G Chromatography
Formulierung	Unconjugated
Formel	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Anwendungsbeschreibung	Flow cytometry: PBMC: Add 10 µl of antibody to 1 million PBMCs in 100 µl PBS. Mix gently and incubate for 15 minutes at 2 to 8°C. Wash twice with PBS, incubate with an appropriate concentration with a fluorochrome conjugated secondary antibody, wash again and analyze or fix with 0.5% v/v of paraformaldehyde in PBS and analyze. WHOLE BLOOD: Add 10 µl of antibody to 100 µl of whole blood. Mix gently and incubate for 15 minutes at room temperature 20°C. Lyse the whole blood. Wash once with PBS, incubate with an appropriate concentration with a fluorochrome conjugated secondary antibody, wash again and analyze or fix with 0.5% v/v of paraformaldehyde in PBS and analyze. See instrument manufacturers instructions for immunofluorescence analysis with a flow cytometer or microscope.