

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Dog Serum NMB-NDS

Artikelname	Dog Serum
Artikelnummer	NMB-NDS
Hersteller Artikelnummer	NDS
Alternativnummer	NMB-NDS
Hersteller	NordicMubio
Kategorie	Zellen/Zellkultur
Spezies Reaktivität	Canine
Produktbeschreibung	As a qualitative serum with an optimized protein content to study the normal protein composition of dog serum. To identify or quantitate a normal serum protein component using a variety of immunodiffusion techniques including immunoelectrophoresis, s...
Puffer	Delipidated, lyophilized, stable reference serum. No preservative added.
Quelle	Pooled normal Dog Serum. The NDS dog serum was obtained from normal, healthy dogs of the Bouvier breed. Regulations are in place in the Netherlands for these dogs to have had some form of vaccination. The Normal Dog Serum (NDS) was derived from dogs that

Anwendungsbeschreibung	Normal serum protein reference sera are shipped lyophilized at ambient temperature. After arrival store at +4C. Prolonged storage may be at -20C. After opening the vial, the lyophilized content is reconstituted by adding 1 ml of sterile distilled water, mixed gently by inversion until complete dissolution is obtained. Allow to stand at ambient temperature for 5-10 minutes to reach equilibrium. Reconstituted serum may be stored frozen.
------------------------	---