

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Goat anti Cat IgA (Fc specific), Clone: [Polyclonal], Monoclonal NMB-GACA/IGA(FC)

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Goat anti Cat IgA (Fc specific), Clone: [Polyclonal], Monoclonal                                                                                                                                                                                                |
| Artikelnummer            | NMB-GACA/IGA(FC)                                                                                                                                                                                                                                                |
| Hersteller Artikelnummer | GACa/IgA(Fc)                                                                                                                                                                                                                                                    |
| Alternativnummer         | NMB-GACA/IGA(FC)                                                                                                                                                                                                                                                |
| Hersteller               | NordicMubio                                                                                                                                                                                                                                                     |
| Wirt                     | Goat                                                                                                                                                                                                                                                            |
| Kategorie                | Antikörper                                                                                                                                                                                                                                                      |
| Produktbeschreibung      | The reactivity of the antiserum is restricted to the Fc part of the IgA molecule. In immunoelectrophoresis and radial immunodiffusion, using various antiserum concentrations against normal cat serum a single precipitin line is obtained which shows a ...   |
| Klonalität               | Monoclonal                                                                                                                                                                                                                                                      |
| Klon-Bezeichnung         | [Polyclonal]                                                                                                                                                                                                                                                    |
| Puffer                   | Delipidated, heat inactivated, lyophilized, stable whole antiserum No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal goat serum. No foreign proteins added. Reconstitute the lyophilized an |
| Quelle                   | Highly purified normal IgA isolated from pooled cat serum. Freund's complete adjuvant is used in the first step of the immunization procedure.                                                                                                                  |

|                        |                                                                                                                                                                                                                                                                                                                        |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formel                 | Delipidated, heat inactivated, lyophilized, stable whole antiserum No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal goat serum. No foreign proteins added.                                                                                        |
| Anwendungsbeschreibung | Precipitation assays. In immunoelectrophoresis use 2 µl serum or equivalent against 120 µl antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 µl antiserum in 3 mm diameter center well and 2 µl serum samples (neat and serially diluted in 2 mm diameter peripheral wells). |