

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Goat anti Guinea Pig IgG2 (subclass specific), conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-GAGP/IGG2/BIO**

|                                  |   |
|----------------------------------|---|
| Artikelname                      | Goat anti Guinea Pig IgG2 (subclass specific), conjugated with Biotin, Clone: [Polyclonal], Monoclonal  |
| Artikelnummer                    | NMB-GAGP/IGG2/BIO   |
| Hersteller Artikelnummer         | GAGp/IgG2/Bio   |
| Alternativnummer                 | NMB-GAGP/IGG2/BIO   |
| Hersteller                       | NordicMubio   |
| Wirt                             | Goat  |
| Kategorie                        | Antikörper  |
| Spezies Reaktivität              | Guinea pig  |
| Konjugation                      | Biotin  |
| Format                           | IgG   |
| Spezifität                       | IgG2 (Fc)   |
| Minimale Kreuzreaktivität (MinX) | no cross-adsorbtion   |
| Produktbeschreibung              | The reactivity of the antiserum is directed to the subclass IgG2. It does not react with IgG1, IgG/Fab fragments, IgM and IgA or any non-Ig protein in guinea pig serum, as tested by immunoelectrophoresis and double radial immunodiffusion. In immunocy... |
| Klonalität                       | Monoclonal  |
| Klon-Bezeichnung                 | [Polyclonal]  |

|                        |   |
|------------------------|---|
| Puffer                 | Biotin-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water, spun do |
| Quelle                 | Purified IgG2 isolated from guinea pig serum. Freund's complete adjuvant is used in the first step of the immunization procedure.   |
| Formulierung           | Biotin  |
| Formel                 | Biotin-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.   |
| Antibody Type          | Secondary Antibody  |
| Anwendungsbeschreibung | ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting.  |