

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Goat anti Swine haptoglobin, Clone: [Polyclonal], Monoclonal NMB-GASW/HP

Artikelname	Goat anti Swine haptoglobin, Clone: [Polyclonal], Monoclonal
Artikelnummer	NMB-GASW/HP
Hersteller Artikelnummer	GASw/Hp
Alternativnummer	NMB-GASW/HP
Hersteller	NordicMubio
Wirt	Goat
Kategorie	Antikörper
Applikation	IP
Spezies Reaktivität	Porcine
Produktbeschreibung	The reactivity of the antiserum is restricted to haptoglobin. In immunoelectrophoresis and radial immunodiffusion (Ouchterlony), using various antiserum concentrations against fresh normal swine plasma a single precipitin line is obtained which shows...
Klonalität	Monoclonal
Klon-Bezeichnung	[Polyclonal]
UniProt	<a href="#">Q8SPS7</a>
Puffer	Delipidated, heat inactivated, lyophilized, stable whole antiserum No preservative added Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal goat serum. No foreign proteins added. Reconstitute the lyophilized ant

Quelle	Haptoglobin is an acute phase protein synthesized in the liver. In immunoelectrophoresis it has an alpha-2 mobility. It forms a group of heterogeneous proteins consisting of four peptide chains, two alpha and two beta chains. The molecular weight is about
Anwendungsbeschreibung	Precipitation assays. In immunoelectrophoresis use 2 ml or equivalent against 120 ml antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 ml in a 3 mm diameter center well and 2 ml serum samples (neat and serially diluted) in 2 mm diameter peripheral wells.