

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Goat anti Swine IgG (heavy and light chains), conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal NMB-GASW/IGG(H+L)/PO**

|                                  |   |
|----------------------------------|---|
| Artikelname                      | Goat anti Swine IgG (heavy and light chains), conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal  |
| Artikelnummer                    | NMB-GASW/IGG(H+L)/PO  |
| Hersteller Artikelnummer         | GASw/IgG(H+L)/PO  |
| Alternativnummer                 | NMB-GASW/IGG(H+L)/PO  |
| Hersteller                       | NordicMubio   |
| Wirt                             | Goat  |
| Kategorie                        | Antikörper  |
| Spezies Reaktivität              | Porcine   |
| Konjugation                      | HRP   |
| Format                           | IgG   |
| Spezifität                       | IgG (H+L)   |
| Minimale Kreuzreaktivität (MinX) | Human, Bovine   |
| Produktbeschreibung              | The reactivity of the antiserum is directed to the Fc and Fab subunits of the IgG molecule. It includes a certain degree of reactivity with other immunoglobulins via the common Fab portion. It does not react with any non-Ig protein in swine serum, as... |
| Klonalität                       | Monoclonal  |
| Klon-Bezeichnung                 | [Polyclonal]  |

|                        |   |
|------------------------|---|
| Puffer                 | Peroxidase-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water, spu |
| Quelle                 | Purified normal IgG isolated from pooled swine serum. Freund's complete adjuvant is used in the first step of the immunization procedure.   |
| Formulierung           | HRP   |
| Formel                 | Peroxidase-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.   |
| Antibody Type          | Secondary Antibody  |
| Anwendungsbeschreibung | ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting.  |