

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Rabbit anti Crotalus durissus terrificus venom Phosphodiesterase I, Clone: [Polyclonal], Monoclonal NMB-NE126/7S**

|                          |   |
|--------------------------|---|
| Artikelname              | Rabbit anti Crotalus durissus terrificus venom Phosphodiesterase I,<br>Clone: [Polyclonal], Monoclonal  |
| Artikelnummer            | NMB-NE126/7S  |
| Hersteller Artikelnummer | NE126/7S  |
| Alternativnummer         | NMB-NE126/7S  |
| Hersteller               | NordicMubio   |
| Wirt                     | Rabbit  |
| Kategorie                | Antikörper  |
| Produktbeschreibung      | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and...   |
| Klonalität               | Monoclonal  |
| Klon-Bezeichnung         | [Polyclonal]  |
| Puffer                   | Purified hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile distilled |
| Quelle                   | Phosphodiesterase I isolated and purified from Crotalus durissus terrificus venom. Freund's complete adjuvant is used in the first step of the immunization procedure.  |

|                        |  |
|------------------------|--|
| Formel                 | Purified hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Anwendungsbeschreibung | Indirect immunofluorescence,ELISA,Dot blot,Immunoblotting.   |