

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Rabbit anti Mouse milk proteins, Clone: [Polyclonal], Monoclonal NMB-RAM/TM

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Rabbit anti Mouse milk proteins, Clone: [Polyclonal], Monoclonal                                                                                                                                                                                                |
| Artikelnummer            | NMB-RAM/TM                                                                                                                                                                                                                                                      |
| Hersteller Artikelnummer | RAM/TM                                                                                                                                                                                                                                                          |
| Alternativnummer         | NMB-RAM/TM                                                                                                                                                                                                                                                      |
| Hersteller               | NordicMubio                                                                                                                                                                                                                                                     |
| Wirt                     | Rabbit                                                                                                                                                                                                                                                          |
| Kategorie                | Antikörper                                                                                                                                                                                                                                                      |
| Produktbeschreibung      | In immunoelectrophoresis against mouse milk precipitation of not less than 7 different proteins may be obtained. However the number of visible precipitin lines varies widely among individual specimens of mouse milk. The total concentration of protein...   |
| Klonalität               | Monoclonal                                                                                                                                                                                                                                                      |
| Klon-Bezeichnung         | [Polyclonal]                                                                                                                                                                                                                                                    |
| Puffer                   | Delipidated, heat inactivated, lyophilized, stable whole antiserum No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal rabbit serum. No foreign proteins added Reconstitute the lyophilized a |
| Quelle                   | Pooled whole mouse milk. Freund's complete adjuvant is used in the first step of the immunization procedure.                                                                                                                                                    |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formel                 | <p>Delipidated, heat inactivated, lyophilized, stable whole antiserum No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal rabbit serum. No foreign proteins added</p>                                                                                                                                                                                                                                                                |
| Anwendungsbeschreibung | <p>Precipitation assays. In immunoelectrophoresis use 2 µl serum or equivalent against 120 µl antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 µl antiserum in 3 mm diameter center well and 2 µl serum samples (neat and serially diluted in 2 mm diameter peripheral wells Different bleedings of the immunized animals are pooled to obtain a broad spectrum balanced against the varying concentrations of the individual serum protein components.</p> |