

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Rabbit anti Photobacterium fischeri NADH-FMN Oxidoreductase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-NE177/BIO**

Artikelname	Rabbit anti Photobacterium fischeri NADH-FMN Oxidoreductase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal
Artikelnummer	NMB-NE177/BIO
Hersteller Artikelnummer	NE177/Bio
Alternativnummer	NMB-NE177/BIO
Hersteller	NordicMubio
Wirt	Rabbit
Kategorie	Antikörper
Konjugation	Biotin
Produktbeschreibung	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and...
Klonalität	Monoclonal
Klon-Bezeichnung	[Polyclonal]
Puffer	Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile dis

Quelle	NADH-FMN oxidoreductase isolated and purified from <i>Photobacterium fischeri</i> . Freund's complete adjuvant is used in the first step of the immunization procedure.
Formulierung	Biotin
Formel	Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Anwendungsbeschreibung	ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting.