

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Anti-Human Eotaxin-2 Antibody, Goat, Polyclonal ABT-ABG10087-U050

|                          |                                                                                                                                                                                                                   |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Anti-Human Eotaxin-2 Antibody, Goat, Polyclonal                                                                                                                                                                   |
| Artikelnummer            | ABT-ABG10087-U050                                                                                                                                                                                                 |
| Hersteller Artikelnummer | ABG10087-U050                                                                                                                                                                                                     |
| Alternativnummer         | ABT-ABG10087-U050-50UG                                                                                                                                                                                            |
| Hersteller               | Abcepta                                                                                                                                                                                                           |
| Wirt                     | Goat                                                                                                                                                                                                              |
| Kategorie                | Antikörper                                                                                                                                                                                                        |
| Applikation              | ELISA, WB                                                                                                                                                                                                         |
| Spezies Reaktivität      | Human                                                                                                                                                                                                             |
| Klonalität               | Polyclonal                                                                                                                                                                                                        |
| Reinheit                 | Produced from sera of goats pre-immunized with highly pure (>98%) recombinant hEotaxin-2. Anti-Human Eotaxin-2 specific antibody was purified by affinity chromatography employing immobilized hEotaxin-2 matrix. |
| Formulierung             | A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.                                                                                                                                            |
| Antibody Type            | Polyclonal Antibody                                                                                                                                                                                               |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>WesternBlot: To detect hEotaxin-2 by Western Blot analysis this antibody can be used at a concentration of 0.1-0.2 µg/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant hEotaxin-2 is 1.5-3.0 ng/lane, under either reducing or non-reducing conditions.. Sandwich: To detect hEotaxin-2 by sandwich ELISA (using 100 µl/well antibody solution) a concentration of 0.5 - 2.0 µg/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with BioGems Biotinylated Anti-Human Eotaxin-2 (60-123BT) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hEotaxin-2.. Reconstitution: Centrifuge vial prior to opening. Reconstitute in sterile water to a concentration of 0.1-1.0 mg/ml.</p> |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|