

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Anti-Murine IFN-gamma Antibody, Rabbit, Polyclonal ABT-ABG10164-U100

|                          |                                                                                                                                                                                                                           |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Anti-Murine IFN-gamma Antibody, Rabbit, Polyclonal                                                                                                                                                                        |
| Artikelnummer            | ABT-ABG10164-U100                                                                                                                                                                                                         |
| Hersteller Artikelnummer | ABG10164-U100                                                                                                                                                                                                             |
| Alternativnummer         | ABT-ABG10164-U100-100UG                                                                                                                                                                                                   |
| Hersteller               | Abcepta                                                                                                                                                                                                                   |
| Wirt                     | Rabbit                                                                                                                                                                                                                    |
| Kategorie                | Antikörper                                                                                                                                                                                                                |
| Applikation              | ELISA, WB                                                                                                                                                                                                                 |
| Spezies Reaktivität      | Mouse                                                                                                                                                                                                                     |
| Klonalität               | Polyclonal                                                                                                                                                                                                                |
| Reinheit                 | Produced from sera of rabbits pre-immunized with highly pure recombinant Murine IFN-gamma. Anti-Murine IFN-gamma specific antibody was purified by affinity chromatography employing immobilized Murine IFN-gamma matrix. |
| Formulierung             | A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.                                                                                                                                                    |
| Antibody Type            | Polyclonal Antibody                                                                                                                                                                                                       |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>WesternBlot: To detect Murine IFN-gamma by Western Blot analysis this antibody can be used at a concentration of 0.1-0.2 µg/ml. When used in conjunction with compatible secondary reagents, the detection limit for recombinant Murine IFN-gamma is 1.5-3.0 ng/lane, under either reducing or non-reducing conditions.. Sandwich: To detect Murine IFN-gamma by sandwich ELISA (using 100 µl/well antibody solution) a concentration of 0.5 - 2.0 µg/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with BioGems Biotinylated Anti-Murine IFN-gamma (61-078BT) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant Murine IFN-gamma.. Neutralization: To yield one-half maximal inhibition [ND&lt;span style=font-size: 16px,&gt;50] of the biological activity of Murine IFN-gamma (1.0 ng/ml), a concentration of 0.05-0.10 ng/ml of this antibody is required.. Reconstitution: Centrifuge vial prior to opening. Reconstitute in sterile water to a concentration of 0.1-1.0 mg/ml.</p> |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|