

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Anti-Human IL-19 Antibody, Rabbit, Polyclonal ABT-ABG10233-U100

|                          |                                                                                                                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Anti-Human IL-19 Antibody, Rabbit, Polyclonal                                                                                                                                                           |
| Artikelnummer            | ABT-ABG10233-U100                                                                                                                                                                                       |
| Hersteller Artikelnummer | ABG10233-U100                                                                                                                                                                                           |
| Alternativnummer         | ABT-ABG10233-U100-100UG                                                                                                                                                                                 |
| Hersteller               | Abcepta                                                                                                                                                                                                 |
| Wirt                     | Rabbit                                                                                                                                                                                                  |
| Kategorie                | Antikörper                                                                                                                                                                                              |
| Applikation              | ELISA, WB                                                                                                                                                                                               |
| Spezies Reaktivität      | Human                                                                                                                                                                                                   |
| Klonalität               | Polyclonal                                                                                                                                                                                              |
| Reinheit                 | Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant hIL-19. Anti-Human IL-19 specific antibody was purified by affinity chromatography employing immobilized hIL-19 matrix. |
| Formulierung             | A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.                                                                                                                                  |
| Antibody Type            | Polyclonal Antibody                                                                                                                                                                                     |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>WesternBlot: To detect hIL-19 by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 µg/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant hIL-19 is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions.. Sandwich: To detect hIL-19 by sandwich ELISA (using 100 µl/well antibody solution) a concentration of 0.5 - 2.0 µg/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with BioGems Biotinylated Anti-Human IL-19 (60-019BT) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hIL-19.. Reconstitution: Centrifuge vial prior to opening. Reconstitute in sterile water to a concentration of 0.1-1.0 mg/ml.</p> |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|