

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **AMACR / p504S (Prostate Cancer Marker) Antibody, IgG2b, Clone: [AMACR/1723], Mouse, Monoclonal NBT-23600-MSM3-P0**

|                          |   |
|--------------------------|---|
| Artikelname              | AMACR / p504S (Prostate Cancer Marker) Antibody, IgG2b, Clone: [AMACR/1723], Mouse, Monoclonal  |
| Artikelnummer            | NBT-23600-MSM3-P0   |
| Hersteller Artikelnummer | 23600-MSM3-P0   |
| Alternativnummer         | NBT-23600-MSM3-P0-20,NBT-23600-MSM3-P0-100  |
| Hersteller               | NeoBiotechnologies  |
| Wirt                     | Mouse   |
| Kategorie                | Antikörper  |
| Applikation              | IHC, WB   |
| Spezies Reaktivität      | Guinea pig, Hamster, Human, Mouse, Rat  |
| Immunogen                | Recombinant full-length human AMACR protein   |
| Produktbeschreibung      | This antibody recognizes a protein of 42kDa, which is identified as AMACR, also known as p504S. It is an enzyme that is involved in bile acid biosynthesis and -oxidation of branched-chain fatty acids. AMACR is essential in lipid metabolism. It is exp... |
| Klonalität               | Monoclonal  |
| Klon-Bezeichnung         | [AMACR/1723]  |
| Molekulargewicht         | 42kDa   |
| Isotyp                   | IgG2b   |

|                        |   |
|------------------------|---|
| NCBI                   | <a href="#">23600</a>   |
| UniProt                | <a href="#">Q9UHK6</a>  |
| Formulierung           | 200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.  |
| Antibody Type          | Monoclonal Antibody   |
| Anwendungsbeschreibung | Western Blot (1-2ug/ml), Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT)(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 min |