

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Cytokeratin 7 (Glandular and Transitional Epithelial Marker) Antibody, IgG1, Clone: [K72.7], Mouse, Monoclonal NBT-3855-MSM2-P0**

|                          |   |
|--------------------------|---|
| Artikelname              | Cytokeratin 7 (Glandular and Transitional Epithelial Marker) Antibody, IgG1, Clone: [K72.7], Mouse, Monoclonal  |
| Artikelnummer            | NBT-3855-MSM2-P0  |
| Hersteller Artikelnummer | 3855-MSM2-P0  |
| Alternativnummer         | NBT-3855-MSM2-P0-20,NBT-3855-MSM2-P0-100  |
| Hersteller               | NeoBiotechnologies  |
| Wirt                     | Mouse   |
| Kategorie                | Antikörper  |
| Applikation              | FC, IF, IHC, WB   |
| Spezies Reaktivität      | Human   |
| Immunogen                | Semi-Purified cytokeratin preparation.  |
| Produktbeschreibung      | It recognizes an intermediate filament protein (IFP) of 55kDa, which is identified as cytokeratin 7. This MAb is highly specific to cytokeratin 7 and shows no cross-reaction with other IFPs. Cytokeratin 7 is a basic cytokeratin, which is found in mos... |
| Klonalität               | Monoclonal  |
| Klon-Bezeichnung         | [K72.7]   |
| Molekulargewicht         | 55kDa   |
| Isotyp                   | IgG1  |

|                        |   |
|------------------------|---|
| NCBI                   | <a href="#">3855</a>  |
| UniProt                | <a href="#">P08729</a>  |
| Formulierung           | 200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.  |
| Antibody Type          | Monoclonal Antibody   |
| Anwendungsbeschreibung | Flow Cytometry (1-2ug/million cells), Immunofluorescence (1-2ug/ml), Western Blot (1-2ug/ml), Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT)(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM E |