

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Mycoplasma genitalium / Mycoplasma hominis PCR kit (Form F) - CE IVD ATB-85-03**

|                          |                                                                                                                                                                                                                                                                    |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Mycoplasma genitalium / Mycoplasma hominis PCR kit (Form F) - CE IVD                                                                                                                                                                                               |
| Artikelnummer            | ATB-85-03                                                                                                                                                                                                                                                          |
| Hersteller Artikelnummer | 85-03                                                                                                                                                                                                                                                              |
| Alternativnummer         | ATB-85-03                                                                                                                                                                                                                                                          |
| Hersteller               | Astra Biotech GmbH                                                                                                                                                                                                                                                 |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                        |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                              |
| Produktbeschreibung      | Qualitative Real-Time PCR detection of Mycoplasma genitalium and Mycoplasma hominis DNA (in tubes)...                                                                                                                                                              |
| Sensitivitaet            | Analytical sensitivity: 1000 copies of M. genitalium, M. hominis genomic DNA per 1 mL                                                                                                                                                                              |
| Proben                   | Urethral swabs, endocervical and vaginal swabs and urine samples                                                                                                                                                                                                   |
| Target-Kategorie         | Mycoplasma genitalium and Mycoplasma hominis                                                                                                                                                                                                                       |
| Anwendungsbeschreibung   | Format F: liquid reaction mix is in 1.5 mL tube). Validated with Qiagen Rotor-Gene Q (also can be used with Corbett Research Rotor-Gene 3000/6000), Bio-Rad CFX96, CFX96 Touch, DNA-Technology DTprime (corresponds to DTlite, SaCycler) and Roche Light Cycler 96 |