

NCBI	2118
Puffer	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Reinheit	Affinity purified
Formulierung	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Application Verdünnung	IF (1:50 - 1:200), IHC (1:50 - 1:200), WB (1:500 - 1:2000)
Anwendungsbeschreibung	The predicted MW is 23kDa/49kDa/53kDa, while the observed MW by Western blot was 55kDa.