

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Caldesmon, HMW (h-Caldesmon) (Smooth Muscle Marker) (rCALD1/820), CF640R conjugate, 0.1mg/mL, Clone: [rCALD1/820], Mouse, Monoclonal BOT-BNC401897-500**

Artikelname	Caldesmon, HMW (h-Caldesmon) (Smooth Muscle Marker) (rCALD1/820), CF640R conjugate, 0.1mg/mL, Clone: [rCALD1/820], Mouse, Monoclonal
Artikelnummer	BOT-BNC401897-500
Hersteller Artikelnummer	BNC401897-500
Alternativnummer	BOT-BNC401897-500-500UL
Hersteller	Biotium
Wirt	Mouse
Kategorie	Antikörper
Applikation	IHC, WB
Spezies Reaktivität	Human
Immunogen	Recombinant full-length human CALD1 protein
Konjugation	CF640R
Produktbeschreibung	Caldesmon HMW is the high molecular weight variant of Caldesmon. Two closely related variants of human caldesmon have been identified which are different in their electrophoretic mobility and cellular distribution. The h-caldesmon variant (120-150 kD...
Klonalität	Monoclonal
Konzentration	0.1 mg/mL

Klon-Bezeichnung	[rCALD1/820]
Molekulargewicht	150 kDa
UniProt	<a href="#">Q05682</a>
Puffer	PBS, 0.1% BSA, 0.05% azide
Quelle	Animal
Anwendungsbeschreibung	Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Does not react with rat Immunohistology (formalin): 0.25-0.5 ug/mL for 30 minutes at RT Staining of formalin-fixed tissues requires boiling tissue sections in 1 mM EDTA pH 7.5-8.5 for 10-20 minutes followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user