

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Villin (GI-Mucosal & Urogenital Brush Border Marker) (rVIL1/1325), CF640R conjugate, 0.1mg/mL, Clone: [rVIL1/1325], Mouse, Monoclonal BOT-BNC402279-500

| | |
|--------------------------|---|
| Artikelname | Villin (GI-Mucosal & Urogenital Brush Border Marker) (rVIL1/1325), CF640R conjugate, 0.1mg/mL, Clone: [rVIL1/1325], Mouse, Monoclonal |
| Artikelnummer | BOT-BNC402279-500 |
| Hersteller Artikelnummer | BNC402279-500 |
| Alternativnummer | BOT-BNC402279-500-500UL |
| Hersteller | Biotium |
| Wirt | Mouse |
| Kategorie | Antikörper |
| Applikation | IHC |
| Spezies Reaktivität | Human |
| Immunogen | Recombinant human Villin fragment of 133 amino acid residues (aa179-311) |
| Konjugation | CF640R |
| Produktbeschreibung | This antibody recognizes a protein of 95 kDa, which is identified as villin. It is a major constituent in the microvilli, which compose the brush border of epithelial cells forming absorptive surfaces of the intestinal and renal proximal tubular epit... |
| Klonalität | Monoclonal |
| Konzentration | 0.1 mg/mL |

| | |
|------------------------|---|
| Klon-Bezeichnung | [rVIL1/1325] |
| Molekulargewicht | 93 kDa |
| UniProt | P09327 |
| Puffer | PBS, 0.1% BSA, 0.05% azide |
| Quelle | Animal |
| Anwendungsbeschreibung | Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Immunohistology (formalin): 1-2 ug/mL for 30 minutes at RT Staining of formalin-fixed tissues requires boiling tissue sections in 10 mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user |