

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Melanoma Marker(HMB45 + M2-7C10 + M2-9E3 + T311), 1mg/mL, Clone: [HMB45 M2-7C10 M2-9E3 T311], Mouse, Monoclonal BOT-BNUM0701-50

Artikelname	Melanoma Marker(HMB45 + M2-7C10 + M2-9E3 + T311), 1mg/mL, Clone: [HMB45 M2-7C10 M2-9E3 T311], Mouse, Monoclonal
Artikelnummer	BOT-BNUM0701-50
Hersteller Artikelnummer	BNUM0701-50
Alternativnummer	BOT-BNUM0701-50-50UL
Hersteller	Biotium
Wirt	Mouse
Kategorie	Antikörper
Applikation	IHC, WB
Spezies Reaktivität	Human
Immunogen	Recombinant hMART-1 protein (M2-7C10, M2-9E3), Recombinant tyrosinase protein (T311), Extract of pigmented melanoma metastases from lymph nodes (HMB45)
Produktbeschreibung	This antibody cocktail recognizes three melanoma-specific proteins, which include MART-1, Tyrosinase and gp100. MART-1 is a newly identified melanocyte differentiation antigen recognized by autologous cytotoxic T lymphocytes. Tyrosinase is one of the...
Klonalität	Monoclonal
Konzentration	1 mg/mL
Klon-Bezeichnung	[HMB45 M2-7C10 M2-9E3 T311]

Molekulargewicht	20-22 kDa (doublet) (MART), 70-80 kDa (Tyrosinase), 90-100 kDa (gp100)
UniProt	Q16655
Puffer	PBS, no BSA, no azide
Quelle	Animal
Anwendungsbeschreibung	Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Immunohistology (Formalin-fixed) 0.5-1.0 ug/mL Staining of formalin-fixed tissues is enhanced by boiling tissue sections in 10 mM citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user