

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Ep-CAM / CD326(Cocktail PAN-EpCAM), 1mg/mL, Clone: [Cocktail PAN-EpCAM], Mouse, Monoclonal BOT-BNUM0969-50

Artikelname	Ep-CAM / CD326(Cocktail PAN-EpCAM), 1mg/mL, Clone: [Cocktail PAN-EpCAM], Mouse, Monoclonal
Artikelnummer	BOT-BNUM0969-50
Hersteller Artikelnummer	BNUM0969-50
Alternativnummer	BOT-BNUM0969-50-50UL
Hersteller	Biotium
Wirt	Mouse
Kategorie	Antikörper
Applikation	Exosome Staining, FC, IF, IHC, WB
Spezies Reaktivität	Human
Immunogen	Recombinant human Ep-CAM protein (full-length and fragments)
Produktbeschreibung	This is a cocktail of four highly specific monoclonal antibodies (EGP40/826, EGP40/837, EGP40/1110, EGP40/1120) that recognize extracellular as well as intracellular domains of the epithelial cellular adhesion molecule (Ep-CAM). Epithelial specific a...
Klonalität	Monoclonal
Konzentration	1 mg/mL
Klon-Bezeichnung	[Cocktail PAN-EpCAM]
Molekulargewicht	40-43 kDa

UniProt	P16422
Puffer	PBS, no BSA, no azide
Quelle	Animal
Anwendungsbeschreibung	Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Immunofluorescence: 1-2 ug/mL Immunohistology formalin-fixed 0.5-1 ug/mL Staining of formalin-fixed tissues requires boiling tissue sections in 10 mM citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes Flow Cytometry 0.5-1 ug/million cells/0.1 mL Western blotting 0.5-1 ug/mL Optimal dilution for a specific application should be determined by user