

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat Anti Pig IgG (H&L)-AbFluor 680 AST-SA1008

Artikelname	Goat Anti Pig IgG (H&L)-AbFluor 680
Artikelnummer	AST-SA1008
Hersteller Artikelnummer	SA1008
Alternativnummer	AST-SA1008-100,AST-SA1008-1
Hersteller	Assay Biotechnology
Wirt	Goat
Kategorie	Antikörper
Applikation	FC, IF, WB
Spezies Reaktivität	Porcine
Konjugation	AbFluor 680
Produktbeschreibung	Assay Biotechnology secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offe...
Quelle	Goat
Reinheit	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Formulierung	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Application Verdünnung	Suggested starting 1:50-1:1000 dilutions for most fluorescent applications, while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR Odyssey equipment would be recommended.
------------------------	---