

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rabbit Anti Guinea Pig IgG (H&L)-AbFluor 680 AST-SA1017

Artikelname	Rabbit Anti Guinea Pig IgG (H&L)-AbFluor 680
Artikelnummer	AST-SA1017
Hersteller Artikelnummer	SA1017
Alternativnummer	AST-SA1017-100,AST-SA1017-1
Hersteller	Assay Biotechnology
Wirt	Rabbit
Kategorie	Antikörper
Applikation	FC, IF, WB
Spezies Reaktivität	Guinea pig
Konjugation	AbFluor 680
Produktbeschreibung	Assay Biotechnology secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offe...
Quelle	Polyclonal, Rabbit IgG
Reinheit	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Formulierung	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Application Verdünnung	Suggested starting 1:50-1:1000 dilutions for most fluorescent applications, while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR Odyssey equipment would be recommended.
------------------------	---