

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Human FABP3 (N-6His), E. coli EKL-PEH0616

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Human FABP3 (N-6His), E. coli                                                                                                                                                                                                                                          |
| Artikelnummer            | EKL-PEH0616                                                                                                                                                                                                                                                                        |
| Hersteller Artikelnummer | PEH0616                                                                                                                                                                                                                                                                            |
| Alternativnummer         | EKL-PEH0616-10UG,EKL-PEH0616-50UG,EKL-PEH0616-500UG                                                                                                                                                                                                                                |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                           |
| Wirt                     | E. coli                                                                                                                                                                                                                                                                            |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                                              |
| Produktbeschreibung      | Fatty Acid Binding Protein 3 (FABP3) is a small cytoplasmic protein (15 kDa) that is released from cardiac myocytes following an ischemic episode. Like the nine other distinct FABPs that have been identified, FABP3 is involved in active fatty acid me...                      |
| Molekulargewicht         | 17.02 KDa                                                                                                                                                                                                                                                                          |
| NCBI                     | <a href="#">05413</a>                                                                                                                                                                                                                                                              |
| Reinheit                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening.Do not mix by vortex or pipetting.It is not recommended to reconstitute to a concentration less than 100µg/ml.Dissolve the lyophilized protein in distilled water.Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |