

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Human PRL2 (C-6His), E. coli EKL-PEH1391

|                          |                                                                                                                                                                                                                                                                                        |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Human PRL2 (C-6His), E. coli                                                                                                                                                                                                                                               |
| Artikelnummer            | EKL-PEH1391                                                                                                                                                                                                                                                                            |
| Hersteller Artikelnummer | PEH1391                                                                                                                                                                                                                                                                                |
| Alternativnummer         | EKL-PEH1391-10UG,EKL-PEH1391-50UG,EKL-PEH1391-500UG                                                                                                                                                                                                                                    |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                               |
| Wirt                     | E. coli                                                                                                                                                                                                                                                                                |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                       |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                                                  |
| Produktbeschreibung      | PTP4A2, also known as PRL2 or PTPCAAX2, is short for Protein tyrosine phosphatase type IVA 2. This protein exists in cell membrane, cytoplasm, endosome and membrane. PTP4A2 is often farnesylated during post-translational modification. Farnesylation is...                         |
| Molekulargewicht         | 20.2 KDa                                                                                                                                                                                                                                                                               |
| NCBI                     | <a href="#">12974</a>                                                                                                                                                                                                                                                                  |
| Reinheit                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                    |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |