

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Human SLAMF6 (C-Fc) EKL-PHH1536

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Human SLAMF6 (C-Fc)                                                                                                                                                                                                                                                    |
| Artikelnummer            | EKL-PHH1536                                                                                                                                                                                                                                                                        |
| Hersteller Artikelnummer | PHH1536                                                                                                                                                                                                                                                                            |
| Alternativnummer         | EKL-PHH1536-10UG,EKL-PHH1536-50UG,EKL-PHH1536-500UG                                                                                                                                                                                                                                |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                           |
| Wirt                     | Human                                                                                                                                                                                                                                                                              |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                                              |
| Produktbeschreibung      | SLAM Family Member 6 (SLAMF6) is a single-pass type I membrane protein that belongs to the SLAM subgroup of the CD2 family. Human SLAMF6/ NTB-A contains a 205 amino acid extracellular domain (ECD) with one Ig-like V-set and one Ig-like C2-set domain,...                      |
| Molekulargewicht         | 49.6 KDa                                                                                                                                                                                                                                                                           |
| NCBI                     | <a href="#">96</a>                                                                                                                                                                                                                                                                 |
| Reinheit                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening.Do not mix by vortex or pipetting.It is not recommended to reconstitute to a concentration less than 100µg/ml.Dissolve the lyophilized protein in distilled water.Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |