

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### GENLISA Human Cluster of Differentiation 82 (CD82) ELISA KBS-KB3011

|                          |                                                                                                                                                                                                                                                              |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | GENLISA Human Cluster of Differentiation 82 (CD82) ELISA                                                                                                                                                                                                     |
| Artikelnummer            | KBS-KB3011                                                                                                                                                                                                                                                   |
| Hersteller Artikelnummer | KB3011                                                                                                                                                                                                                                                       |
| Alternativnummer         | KBS-KB3011-1                                                                                                                                                                                                                                                 |
| Hersteller               | Krishgen Biosystems                                                                                                                                                                                                                                          |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                  |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                        |
| Produktbeschreibung      | The GENLISA Human CD82 ELISA includes features like: - Ready to use protocol with break-apart wells for ease of use - Standardisation and High Reproducibility - Lot to Lot Consistency - Accuracy and Precision Validated against seven points for a "GO... |
| Detektionsbereich        | 78.125 - 5000 pg/ml                                                                                                                                                                                                                                          |
| Sensitivitaet            | 70 pg/ml                                                                                                                                                                                                                                                     |
| Proben                   | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                           |
| Anwendungsbeschreibung   | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                   |