

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### GENLISA Human Dc Specific Intercellular Adhesion Molecule 3-Grabbing Nonintegrin (DC-SIGN / CD209) ELISA KBS-KB3089

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | GENLISA Human Dc Specific Intercellular Adhesion Molecule 3-Grabbing Nonintegrin (DC-SIGN / CD209) ELISA                                                                                                                                                      |
| Artikelnummer            | KBS-KB3089                                                                                                                                                                                                                                                    |
| Hersteller Artikelnummer | KB3089                                                                                                                                                                                                                                                        |
| Alternativnummer         | KBS-KB3089-1                                                                                                                                                                                                                                                  |
| Hersteller               | Krishgen Biosystems                                                                                                                                                                                                                                           |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                         |
| Produktbeschreibung      | The GENLISA Human Dc Specific Intercellular Adhesion Molecule 3-Grabbing Nonintegrin (DC-SIGN / CD209) ELISA includes features like: - Ready to use protocol with break-apart wells for ease of use - Standardisation and High Reproducibility - Lot to Lo... |
| Detektionsbereich        | 23.44 - 1500 pg/ml                                                                                                                                                                                                                                            |
| Sensitivitaet            | 18 pg/ml                                                                                                                                                                                                                                                      |
| Proben                   | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                            |
| Anwendungsbeschreibung   | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                    |