

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

BNP and NT-proBNP Free Plasma Preis auf Anfrage RAY-MD-17-0049P

Artikelname	BNP and NT-proBNP Free Plasma Preis auf Anfrage
Artikelnummer	RAY-MD-17-0049P
Hersteller Artikelnummer	MD-17-0049P
Alternativnummer	RAY-MD-17-0049P
Hersteller	RayBiotech
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Produktbeschreibung	BNP and NT-proBNP Free Plasma...
Konzentration	Not Applicable
Expression System	Human Plasma
Reinheit	Immunoaffinity chromatography. Gel for affinity column was prepared by coupling anti-BNP and anti-NTproBNP monoclonal antibodies with BrCN-activated Sepharose CL 4B.
Formulierung	Affinity Purified, Liquid
Formel	Not Applicable
Application Verdünnung	CON Can be used for preparation of any proBNP-related standards and calibrators. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not nec