

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **PerCP/Cyanine5.5 Anti-Human CD371 Antibody[50C1], PerCP/Cy5.5, Mouse, Monoclonal ELA-AN00320J**

|                          |                                                                                  |
|--------------------------|----------------------------------------------------------------------------------|
| Artikelname              | PerCP/Cyanine5.5 Anti-Human CD371 Antibody[50C1], PerCP/Cy5.5, Mouse, Monoclonal |
| Artikelnummer            | ELA-AN00320J                                                                     |
| Hersteller Artikelnummer | AN00320J                                                                         |
| Alternativnummer         | ELA-AN00320J-100TESTS,ELA-AN00320J-100TESTSX2,ELA-AN00320J-20TESTS               |
| Hersteller               | Elabscience                                                                      |
| Wirt                     | Mouse                                                                            |
| Kategorie                | Antikörper                                                                       |
| Applikation              | FC                                                                               |
| Spezies Reaktivität      | Human                                                                            |
| Konjugation              | PerCP/Cy5.5                                                                      |
| Klonalität               | Monoclonal                                                                       |
| Konzentration            | 5 µL/Test                                                                        |
| Klon-Bezeichnung         | [50C1]                                                                           |
| NCBI                     | <a href="#">160364</a>                                                           |
| UniProt                  | <a href="#">Q5QGZ9</a>                                                           |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                           |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Puffer                 | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.                                                                                                                                                                                                                                                                                                               |
| Target-Kategorie       | CD371                                                                                                                                                                                                                                                                                                                                                                                                     |
| Anwendungsbeschreibung | Each lot of this antibody is quality control tested by flow cytometric analysis. <b style="color:0070ca;">The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.</b> |