

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rabbit CXCL8 / IL8 protein, His tag (active), Unconjugated GTX00339-PRO

Artikelname	Rabbit CXCL8 / IL8 protein, His tag (active), Unconjugated
Artikelnummer	GTX00339-PRO
Hersteller Artikelnummer	GTX00339-pro
Alternativnummer	GTX00339-PRO-10
Hersteller	GeneTex
Kategorie	Proteine/Peptide
Applikation	FA
Spezies Reaktivität	Rabbit
Konjugation	Unconjugated
NCBI	100009129
UniProt	P19874
Puffer	Reconstitute with 20mM Tris and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose, ProClin 300.
Expression System	E. coli
Formulierung	Lyophilized powder
Sequenz	N-terminal His-Tag, Ala23~Ser101 (NP_001075762.1)

Anwendungsbeschreibung

Interleukin 8 (IL8 or chemokine (C-X-C motif) ligand 8, CXCL8) is a chemokine produced by macrophages and other cell types such as epithelial cells, airway smooth muscle cells and endothelial cells. Thus, chemotaxis assay used 24-well microchemotaxis system was undertaken to detect the chemotactic effect of IL8 on the human T-lymphocyte leukemia cell line Jurkat. Briefly, Jurkat cells were seeded into the upper chambers (100 μ l cell suspension, 10^6 cells/ml in RPMI-1640 with FBS free) and IL8 (0.1 ng/ml, 1 ng/ml and 10 ng/ml diluted separately in serum free RPMI-1640) was added in lower chamber with a polycarbonate filter (8 μ m pore size) used to separate the two compartments. After incubation at 37°C with 5% CO₂ for 1h, the filter was removed, then cells in low chamber were observed by inverted microscope at low magnification ($\times 100$) and the number of migrated cells were counted at high magnification ($\times 400$) randomly (five fields for each filter). The migrated Jurkat cells in low chamber at low magnification ($\times 100$) randomly, and the migrated cells were counted at high magnification ($\times 400$). IL8 is able to induce migration of Jurkat cells, and the optimum chemotaxis of IL8 occurs at 1~10 ng/ml.