

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### DARC blocking peptide GTX89467-PEP

|                          |                                                                                                                                              |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | DARC blocking peptide                                                                                                                        |
| Artikelnummer            | GTX89467-PEP                                                                                                                                 |
| Hersteller Artikelnummer | GTX89467-PEP                                                                                                                                 |
| Alternativnummer         | GTX89467-PEP-100                                                                                                                             |
| Hersteller               | GeneTex                                                                                                                                      |
| Kategorie                | Proteine/Peptide                                                                                                                             |
| Spezies Reaktivität      | Human                                                                                                                                        |
| Sensitivitaet            | The peptide was used in the production of GTX89467.                                                                                          |
| NCBI                     | <a href="#">2532</a>                                                                                                                         |
| Puffer                   | Reconstitute with 200µl distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H2O, no preservatives. |
| Formulierung             | Lyophilized powder                                                                                                                           |
| Sequenz                  | HRAELSPSTENSSQLDFED-C                                                                                                                        |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>The antibodies(GTX89467) come at 0.5mg/ml and peptide as a 100ug pellet. When peptides are reconstituted in 200ul water, the concentration would be also 0.5mg/ml. To start, the best ratio would be 1:1 (which means molar excess of peptides relative to antibodies when identical volumes are mixed).&lt;br/&gt;For Peptide Blocking:&lt;br/&gt;Please mix equal volumes of peptide and antibody at the required dilution and leave at ambient temperature. Best is to have two identical blots, to be incubated with equal amount of antibodies, but one with the antibodies pre-adsorbed to the peptide for 20min. Then incubate and develop in parallel the two blots.</p> |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|