

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Cytochrome C Antibody, Clone: [SPM389], Mouse, Monoclonal NSJ-V2787IHC-7ML

|                          |   |
|--------------------------|---|
| Artikelname              | Cytochrome C Antibody, Clone: [SPM389], Mouse, Monoclonal   |
| Artikelnummer            | NSJ-V2787IHC-7ML  |
| Hersteller Artikelnummer | V2787IHC-7ML  |
| Alternativnummer         | NSJ-V2787IHC-7ML  |
| Hersteller               | NSJ Bioreagents   |
| Wirt                     | Mouse   |
| Kategorie                | Antikörper  |
| Applikation              | IHC-P   |
| Spezies Reaktivität      | Human   |
| Immunogen                | Amino acids 66-104 of pigeon Cytochrome C were used as the immunogen for the Cytochrome C antibody. It recognizes an epitope within amino acids 93-104 of pigeon Cytochrome C.  |
| Produktbeschreibung      | Cytochrome C is a well-characterized mobile electron transport protein that is essential to energy conversion in all aerobic organisms. In mammalian cells, this highly conserved protein is normally localized to the mitochondrial inter-membrane space.... |
| Klonalität               | Monoclonal  |
| Klon-Bezeichnung         | [SPM389]  |
| UniProt                  | <a href="#">P99999</a>  |
| Reinheit                 | Protein G affinity chromatography   |

|                        |  |
|------------------------|--|
| Formulierung           | Prediluted in 1X PBS, 0.1 mg/ml BSA (US sourced), 0.05% sodium azide, *For IHC use only*   |
| Antibody Type          | Primary Antibody   |
| Application Verdünnung | The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.  |
| Anwendungsbeschreibung | Optimal dilution of the Cytochrome C antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in pH 9 10mM Tris with 1mM EDTA for 10-20 min followed by cooling at RT for 20 min2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min. |