

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Melanoma Marker Antibody Cocktail (MART-1 & Tyrosinase & gp100), Clone: [DT101+ BC199 + T311 + HMB45], Mouse, Monoclonal NSJ-V8328-20UG

Artikelname	Melanoma Marker Antibody Cocktail (MART-1 & Tyrosinase & gp100), Clone: [DT101+ BC199 + T311 + HMB45], Mouse, Monoclonal
Artikelnummer	NSJ-V8328-20UG
Hersteller Artikelnummer	V8328-20UG
Alternativnummer	NSJ-V8328-20UG
Hersteller	NSJ Bioreagents
Wirt	Mouse
Kategorie	Antikörper
Applikation	IHC-P
Spezies Reaktivität	Human
Immunogen	Recombinant hMART-1 protein (DT101 & BC199), extract of pigmented melanoma metastases from lymph nodes (HMB45) and recombinant tyrosinase protein (T311) were used as the immunogen for this Melanoma Marker antibody cocktail.
Produktbeschreibung	This antibody cocktail recognizes three melanoma-specific proteins, which include MART-1, Tyrosinase and gp100. MART-1 is a newly identified melanocyte differentiation antigen recognized by autologous cytotoxic T lymphocytes. Tyrosinase is one of the...
Klonalität	Monoclonal
Klon-Bezeichnung	[DT101+ BC199 + T311 + HMB45]

Reinheit	Protein G affinity chromatography
Formulierung	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Antibody Type	Primary Antibody
Application Verdünnung	Immunohistochemistry (FFPE): 1-2ug/ml
Anwendungsbeschreibung	Optimal dilution of the Melanoma Marker antibody cocktail should be determined by the researcher.