

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rat Trypsinogen activation peptide (TAP) ELISA Kit ASC-KTE100184

Artikelname	Rat Trypsinogen activation peptide (TAP) ELISA Kit
Artikelnummer	ASC-KTE100184
Hersteller Artikelnummer	KTE100184
Alternativnummer	ASC-KTE100184-48, ASC-KTE100184-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Rat
Produktbeschreibung	This Rat Trypsinogen activation peptide (TAP) ELISA Kit employs a two-site sandwich ELISA to quantitate TAP....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	TAP
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Rat Trypsinogen activation peptide (TAP) ELISA Kit employs a two-site sandwich ELISA to quantitate TAP in samples. An antibody specific for TAP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TAP present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TAP is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TAP bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	--