

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Rat Amyloid Beta Peptide 1-42 (Abeta 1-42) ELISA Kit ASC-KTE101180

|                          |                                                                                |
|--------------------------|--------------------------------------------------------------------------------|
| Artikelname              | Rat Amyloid Beta Peptide 1-42 (Abeta 1-42) ELISA Kit                           |
| Artikelnummer            | ASC-KTE101180                                                                  |
| Hersteller Artikelnummer | KTE101180                                                                      |
| Alternativnummer         | ASC-KTE101180-48, ASC-KTE101180-96                                             |
| Hersteller               | Abbkine Scientific                                                             |
| Kategorie                | Kits/Assays                                                                    |
| Applikation              | ELISA                                                                          |
| Spezies Reaktivität      | Rat                                                                            |
| Produktbeschreibung      | Null employs a two-site sandwich ELISA to quantitate Abeta 1-42 in samples.... |
| Detektionsbereich        | 3 µg/L-48 µg/L                                                                 |
| Sensitivitaet            | 0.1µg/L                                                                        |
| Tag                      | Abeta 1-42                                                                     |
| Proben                   | Cell culture supernatantsSerumPlasmaOther biological fluids                    |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>Rat Amyloid Beta Peptide 1-42 (Abeta 1-42) ELISA Kit employs a two-site sandwich ELISA to quantitate Abeta 1-42 in samples. An antibody specific for Abeta 1-42 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any Abeta 1-42 present is bound by the immobilized antibody. After removing any unbound substances, HRP-Conjugate Rat Abeta 1-42 detection antibody is added to the wells. Following a wash to remove any unbound HRP reagent, a Chromogen solution is added to the wells and color develops in proportion to the amount of Abeta 1-42 bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|