

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Human Sodium- and chloride-dependent taurine transporter (SLC6A6/TAUT) ELISA Kit ASC-KTE60592

Artikelname	Human Sodium- and chloride-dependent taurine transporter (SLC6A6/TAUT) ELISA Kit
Artikelnummer	ASC-KTE60592
Hersteller Artikelnummer	KTE60592
Alternativnummer	ASC-KTE60592-48, ASC-KTE60592-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human Sodium- and chloride-dependent taurine transporter (SLC6A6/TAUT) ELISA Kit employs a two-site sandwich ELISA to quantitate SLC6A6....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	SLC6A6
NCBI	6533
UniProt	P31641
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Human Sodium- and chloride-dependent taurine transporter (SLC6A6/TAUT) ELISA Kit employs a two-site sandwich ELISA to quantitate SLC6A6 in samples. An antibody specific for SLC6A6 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any SLC6A6 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for SLC6A6 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of SLC6A6 bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	---