

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Human Putative 2-oxo-4-hydroxy-4-carboxy-5-ureidoimidazoline decarboxylase (PRHOXNB) ELISA Kit ASC-KTE61163

Artikelname	Human Putative 2-oxo-4-hydroxy-4-carboxy-5-ureidoimidazoline decarboxylase (PRHOXNB) ELISA Kit
Artikelnummer	ASC-KTE61163
Hersteller Artikelnummer	KTE61163
Alternativnummer	ASC-KTE61163-48, ASC-KTE61163-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human Putative 2-oxo-4-hydroxy-4-carboxy-5-ureidoimidazoline decarboxylase (PRHOXNB) ELISA Kit employs a two-site sandwich ELISA to quantitate PRHOXNB....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	PRHOXNB
NCBI	646625
UniProt	A6NGE7
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Human Putative 2-oxo-4-hydroxy-4-carboxy-5-ureidoimidazoline decarboxylase (PRHOXNB) ELISA Kit employs a two-site sandwich ELISA to quantitate PRHOXNB in samples. An antibody specific for PRHOXNB has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any PRHOXNB present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for PRHOXNB is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of PRHOXNB bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	--