

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Human Oxytocin (OT) ELISA Kit ASC-KTE62212

|                          |                                                                                                       |
|--------------------------|-------------------------------------------------------------------------------------------------------|
| Artikelname              | Human Oxytocin (OT) ELISA Kit                                                                         |
| Artikelnummer            | ASC-KTE62212                                                                                          |
| Hersteller Artikelnummer | KTE62212                                                                                              |
| Alternativnummer         | ASC-KTE62212-48, ASC-KTE62212-96                                                                      |
| Hersteller               | Abbkine Scientific                                                                                    |
| Kategorie                | Kits/Assays                                                                                           |
| Applikation              | ELISA                                                                                                 |
| Spezies Reaktivität      | Human                                                                                                 |
| Produktbeschreibung      | Human Oxytocin (OT) ELISA Kit employs a two-site sandwich ELISA to quantitate Oxytocin in samples.... |
| Detektionsbereich        | 2.5 ng/L-40 ng/L                                                                                      |
| Sensitivitaet            | 0.1 ng/L                                                                                              |
| Tag                      | Oxytocin                                                                                              |
| Proben                   | Cell culture supernatantsSerumPlasmaOther biological fluids                                           |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>Human Oxytocin (OT) ELISA Kit employs a two-site sandwich ELISA to quantitate Oxytocin in samples. An antibody specific for Oxytocin has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any Oxytocin present is bound by the immobilized antibody. After removing any unbound substances, HRP-Conjugate Human Oxytocin detection antibody is added to the wells. Following a wash to remove any unbound HRP reagent, a Chromogen solution is added to the wells and color develops in proportion to the amount of Oxytocin bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|