

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse Vasoactive Intestinal Peptide (VIP) ELISA Kit ASC-KTE70020

Artikelname	Mouse Vasoactive Intestinal Peptide (VIP) ELISA Kit
Artikelnummer	ASC-KTE70020
Hersteller Artikelnummer	KTE70020
Alternativnummer	ASC-KTE70020-48, ASC-KTE70020-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Mouse
Produktbeschreibung	This Mouse Vasoactive Intestinal Peptide (VIP) ELISA Kit employs a two-site sandwich ELISA to quantitate VIP....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	VIP
NCBI	22353
UniProt	P32648
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Mouse Vasoactive Intestinal Peptide (VIP) ELISA Kit employs a two-site sandwich ELISA to quantitate VIP in samples. An antibody specific for VIP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyVIP present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for VIP is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of VIP bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	--