

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse Toll-interacting protein (TOLLIP) ELISA Kit ASC-KTE70080

Artikelname	Mouse Toll-interacting protein (TOLLIP) ELISA Kit
Artikelnummer	ASC-KTE70080
Hersteller Artikelnummer	KTE70080
Alternativnummer	ASC-KTE70080-48, ASC-KTE70080-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Mouse
Produktbeschreibung	This Mouse Toll-interacting protein (TOLLIP) ELISA Kit employs a two-site sandwich ELISA to quantitate TOLLIP....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	TOLLIP
NCBI	54473
UniProt	Q9QZ06
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Mouse Toll-interacting protein (TOLLIP) ELISA Kit employs a two-site sandwich ELISA to quantitate TOLLIP in samples. An antibody specific for TOLLIP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TOLLIP present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TOLLIP is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TOLLIP bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	--