

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse Tenascin-N (TNN) ELISA Kit ASC-KTE70103

Artikelname	Mouse Tenascin-N (TNN) ELISA Kit
Artikelnummer	ASC-KTE70103
Hersteller Artikelnummer	KTE70103
Alternativnummer	ASC-KTE70103-48, ASC-KTE70103-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Mouse
Produktbeschreibung	This Mouse Tenascin-N (TNN) ELISA Kit employs a two-site sandwich ELISA to quantitate TNN....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	TNN
NCBI	329278
UniProt	Q80Z71
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Mouse Tenascin-N (TNN) ELISA Kit employs a two-site sandwich ELISA to quantitate TNN in samples. An antibody specific for TNN has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TNN present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNN is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNN bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	--