

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Mouse Tumor necrosis factor alpha-induced protein 8 (TNFAIP8) ELISA Kit ASC-KTE70134

|                          |                                                                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Mouse Tumor necrosis factor alpha-induced protein 8 (TNFAIP8) ELISA Kit                                                                  |
| Artikelnummer            | ASC-KTE70134                                                                                                                             |
| Hersteller Artikelnummer | KTE70134                                                                                                                                 |
| Alternativnummer         | ASC-KTE70134-48, ASC-KTE70134-96                                                                                                         |
| Hersteller               | Abbkine Scientific                                                                                                                       |
| Kategorie                | Kits/Assays                                                                                                                              |
| Applikation              | ELISA                                                                                                                                    |
| Spezies Reaktivität      | Mouse                                                                                                                                    |
| Produktbeschreibung      | This Mouse Tumor necrosis factor alpha-induced protein 8 (TNFAIP8) ELISA Kit employs a two-site sandwich ELISA to quantitate TNFAIP8.... |
| Detektionsbereich        | Please inquire                                                                                                                           |
| Sensitivitaet            | Please inquire                                                                                                                           |
| Tag                      | TNFAIP8                                                                                                                                  |
| NCBI                     | <a href="#">106869</a>                                                                                                                   |
| UniProt                  | <a href="#">Q921Z5</a>                                                                                                                   |
| Proben                   | Cell culture supernatantsSerumPlasmaOther biological fluids                                                                              |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>This Mouse Tumor necrosis factor alpha-induced protein 8 (TNFAIP8) ELISA Kit employs a two-site sandwich ELISA to quantitate TNFAIP8 in samples. An antibody specific for TNFAIP8 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyTNFAIP8 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNFAIP8 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNFAIP8 bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|