

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse TRAF family member-associated NF-kappa-B activator (TANK) ELISA Kit ASC-KTE70335

Artikelname	Mouse TRAF family member-associated NF-kappa-B activator (TANK) ELISA Kit
Artikelnummer	ASC-KTE70335
Hersteller Artikelnummer	KTE70335
Alternativnummer	ASC-KTE70335-48, ASC-KTE70335-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Mouse
Produktbeschreibung	This Mouse TRAF family member-associated NF-kappa-B activator (TANK) ELISA Kit employs a two-site sandwich ELISA to quantitate TANK....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	TANK
NCBI	21353
UniProt	P70347
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Mouse TRAF family member-associated NF-kappa-B activator (TANK) ELISA Kit employs a two-site sandwich ELISA to quantitate TANK in samples. An antibody specific for TANK has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TANK present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TANK is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TANK bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	--