

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Severe acute respiratory syndrome coronavirus 2 Spike glycoprotein Protein BYT-ORB1477212

Artikelname	Severe acute respiratory syndrome coronavirus 2 Spike glycoprotein Protein
Artikelnummer	BYT-ORB1477212
Hersteller Artikelnummer	orb1477212
Alternativnummer	BYT-ORB1477212-1,BYT-ORB1477212-100,BYT-ORB1477212-20
Hersteller	Biorbyt
Kategorie	Proteine/Peptide
Produktbeschreibung	This Severe acute respiratory syndrome coronavirus 2 Spike glycoprotein Protein spans the amino acid sequence from region 319-541aa. Purity: Greater than 90% as determined by SDS-PAGE....
Molekulargewicht	28.3 kDa
UniProt	P0DTC2
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)
Reinheit	Greater than 90% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder

Sequenz	RVQPTESIVRFPNITNLCPFGEVFNATRFASVYAWNRRKRISNCVADYSVLYNSA SFSTFKCYGVSPTKLNDLCFTNVYADSFVIRGDEVQRQIAPGQTGKIADYNYKLP DDFTGCVIAWNSNNLDSKVGGNYNYLYRLFRKSNLKPFERDISTEIYQAGSTP CNGVEGFNCYFPLQSYGFQPTNGVGYQPYRVVVLSEFELLHAPATVCGPKKST NLVKNKCVNF
Anwendungsbeschreibung	Biological Origin: Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2). Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference