

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Reovirus type 3 Outer capsid protein sigma-3 Protein BYT-ORB1477218

Artikelname	Reovirus type 3 Outer capsid protein sigma-3 Protein
Artikelnummer	BYT-ORB1477218
Hersteller Artikelnummer	orb1477218
Alternativnummer	BYT-ORB1477218-1,BYT-ORB1477218-100,BYT-ORB1477218-20
Hersteller	Biorbyt
Kategorie	Proteine/Peptide
Produktbeschreibung	This Reovirus type 3 Outer capsid protein sigma-3 Protein spans the amino acid sequence from region 1-365aa. Purity: Greater than 90% as determined by SDS-PAGE....
Molekulargewicht	46.7 kDa
UniProt	<a href="#">P03527</a>
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	Reovirus type 3 (strain Dearing) (T3D) (Mammalian orthoreovirus 3)
Reinheit	Greater than 90% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder

Sequenz	MEVCLPNGHQVVDLINNAFEGRVSIYSAQEGWDKTISAQPDMMVCGGAVVC MHCLGVVGS LQRK LKHLPHHRCNQQRHQDYVDVQFADRVT AHWKRGMLS FVAQMHEMMNDVSPDDLDRVRTEGGSLVELNWLQVDPNSMFRSIHSSWTD PLQVVDDLDTKLDQYWTALNLMIDSSDLIPNFMMRDPSHAFNGVKLGGDAR QTQFSRTFDSRSSLEWGMVYDYSELEHDP SKGRAYRKELVTPARDFGHFGL S
Anwendungsbeschreibung	Biological Origin: Reovirus type 3 (strain Dearing) (T3D) (Mammalian orthoreovirus 3). Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference