

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rat Excitatory Amino Acid Transporter 3 (EAAT3) ELISA Kit BYT-ORB1947837

Artikelname	Rat Excitatory Amino Acid Transporter 3 (EAAT3) ELISA Kit
Artikelnummer	BYT-ORB1947837
Hersteller Artikelnummer	orb1947837
Alternativnummer	BYT-ORB1947837-48,BYT-ORB1947837-96
Hersteller	Biorbyt
Kategorie	Kits/Assays
Spezies Reaktivität	Rat
Produktbeschreibung	Rat Excitatory Amino Acid Transporter 3 (EAAT3) ELISA Kit...
Detektionsbereich	0.63-40ng/mL
Sensitivitaet	0.38 ng/mL
Proben	serum, plasma, Tissue homogenate and Other biological samples
Target-Kategorie	EAAT3

Anwendungsbeschreibung

Application Notes: This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Rat EAAT3. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Rat EAAT3 and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Rat EAAT3, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm. The OD value is proportional to the concentration of Rat EAAT3. You can calculate the concentration of Rat EAAT3 in the samples by comparing the OD of the samples to the standard curve