

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Severe acute respiratory syndrome coronavirus 2 Replicase polyprotein 1ab (rep) (P4714L), partial BYT-ORB1096040

Artikelname	Recombinant Severe acute respiratory syndrome coronavirus 2 Replicase polyprotein 1ab (rep) (P4714L), partial
Artikelnummer	BYT-ORB1096040
Hersteller Artikelnummer	orb1096040
Alternativnummer	BYT-ORB1096040-20,BYT-ORB1096040-100,BYT-ORB1096040-1
Hersteller	Biorbyt
Kategorie	Proteine/Peptide
Produktbeschreibung	This Recombinant Severe acute respiratory syndrome coronavirus 2 RNA-dependent RNA polymerase (NSP12) spans the amino acid sequence from region 1-932aa. Purity: Greater than 85% as determined by SDS-PAGE....
Molekulargewicht	108.2 kDa
UniProt	<a href="#">P0DTD1</a>
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder

Sequenz	SADAQSFLNRVCGVSAARLTPCGTGTSTDVVYRAFDIYNDKVAGFAKFLKTN CCRFQEKDEDDNLIDSYFVVKRHTFSNYQHEETIYNLLKDCPAVAKHDFFKFRI DGDMVPHISRQRLTKYTMADLVYALRHFDEGNCDTLKEILVTYNCCDDDDYFN KKDWYDFVENPDILRVYANLGERVRQALLKTVQFCDAMRNAGIVGVLTLDNQ DLNGNWYDFGDFIQTTTPGSGVPVVDSSYSSLLMPILTLTRALTAES
Anwendungsbeschreibung	Biological Origin: Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2). Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference