

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant *Macrococcus caseolyticus* Oleate hydratase (MCCL_0076), partial BYT-ORB1096062

| | |
|--------------------------|--|
| Artikelname | Recombinant <i>Macrococcus caseolyticus</i> Oleate hydratase (MCCL_0076), partial |
| Artikelnummer | BYT-ORB1096062 |
| Hersteller Artikelnummer | orb1096062 |
| Alternativnummer | BYT-ORB1096062-20,BYT-ORB1096062-100,BYT-ORB1096062-1 |
| Hersteller | Biorbyt |
| Kategorie | Proteine/Peptide |
| Produktbeschreibung | This Recombinant <i>Macrococcus caseolyticus</i> Oleate hydratase (MCCL_0076), partial spans the amino acid sequence from region 1-367aa. Purity: Greater than 85% as determined by SDS-PAGE.... |
| Molekulargewicht | 45.5 kDa |
| UniProt | B9E972 |
| Puffer | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0. |
| Quelle | <i>Macrococcus caseolyticus</i> (strain JCSC5402) |
| Reinheit | Greater than 85% as determined by SDS-PAGE. |
| Formulierung | Liquid or Lyophilized powder |

| | |
|------------------------|--|
| Sequenz | MYYSNGNYEAFARPKKPEGVDNKSAYLVGSGLASLAAASFLIRDGQMKGENI HILEELDLPGGSLDGILNPERGYIMRGGREMENHFECLEWDLFRSVPSLEVEDA SVLDEFYWLNKEDPNYSKCRVIENRGQRLESDGKMTLTKKANKEIIQLCLMKE EQLNDVKISDVFSKDFLDSNFWIYWKTMFAPWHSAMEMRRYLMRFIHHIG GLADFSALKFTKFNQFESLVMPLEHLKAKNVTFEYGVTVKNIQV |
| Anwendungsbeschreibung | <p>Biological Origin: <i>Macrocooccus caseolyticus</i> (strain JCSC5402).</p> <p>Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference</p> |