

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Reovirus type 3 Outer capsid protein lambda-2 (L2), partial BYT-ORB1096393

Artikelname	Recombinant Reovirus type 3 Outer capsid protein lambda-2 (L2), partial
Artikelnummer	BYT-ORB1096393
Hersteller Artikelnummer	orb1096393
Alternativnummer	BYT-ORB1096393-20,BYT-ORB1096393-100,BYT-ORB1096393-1
Hersteller	Biorbyt
Kategorie	Proteine/Peptide
Produktbeschreibung	This Recombinant Reovirus type 3 Outer capsid protein lambda-2 (L2), partial spans the amino acid sequence from region 1-200aa. Purity: Greater than 90% as determined by SDS-PAGE....
Molekulargewicht	30.2 kDa
UniProt	P11079
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	Reovirus type 3 (strain Dearing) (T3D) (Mammalian orthoreovirus 3)
Reinheit	Greater than 90% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder

Sequenz	MANVWGVRLADSLSSPTIETRTRQYTLHDLCSLDANPGREPWKPLRNQRT NNIVAVQLFRPLQGLVLDTQLYGFPGAFDDWERFMREKLRVLKYEVLRIYPISN YSNEHVNVFVANALVGAFLSNQAFYDLLPLLIINDTMIGDLLGTGASLSQFFQS HGDVLEVAAGRKYLQMENYSNDDDDPPLFAKDLSDYAKAFY
Anwendungsbeschreibung	Biological Origin: Reovirus type 3 (strain Dearing) (T3D) (Mammalian orthoreovirus 3). Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference