

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Thermus aquaticus DNA mismatch repair protein mutL (mutL), partial BYT-ORB1096441

| | |
|--------------------------|--|
| Artikelname | Recombinant Thermus aquaticus DNA mismatch repair protein mutL (mutL), partial |
| Artikelnummer | BYT-ORB1096441 |
| Hersteller Artikelnummer | orb1096441 |
| Alternativnummer | BYT-ORB1096441-20,BYT-ORB1096441-100,BYT-ORB1096441-1 |
| Hersteller | Biorbyt |
| Kategorie | Proteine/Peptide |
| Produktbeschreibung | This Recombinant Thermus aquaticus DNA mismatch repair protein mutL (mutL), partial spans the amino acid sequence from region 186-312aa. Purity: Greater than 90% as determined by SDS-PAGE.... |
| Molekulargewicht | 21.3 kDa |
| UniProt | P96082 |
| Puffer | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0. |
| Quelle | Thermus aquaticus |
| Reinheit | Greater than 90% as determined by SDS-PAGE. |
| Formulierung | Liquid or Lyophilized powder |

| | |
|------------------------|---|
| Sequenz | LFPGAGLKEAARQAFGGLLAERLFPLEKGGAFLEGLLTGPQVSRTRPDLLFLA VNGRPVALPEGVLRVRRRAYRELLPEGHYPVGVNLNLSLPPGAYRLRLDARKEE VALSKEAEAFLEEALIEEFR |
| Anwendungsbeschreibung | Biological Origin: <i>Thermus aquaticus</i> . Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference |